



RhD typing

Bed side test

Antibody tests

Practice for IV year medical students

Zita Csernus MD

National Blood Transfusion Service

Blood Transfusion Centre Pécs




www.ovsz.hu

Rh Blood Group System

Discovery: 1940 K. Landsteiner and A. Wiener

Experience: Rabbits were injected with the red blood cells of the Rhesus monkey. The injection caused an antigenic reaction in the serum component of rabbit blood.



When blood from **humans was tested** with the rabbit serum, the red blood cells of 85% of the humans tested agglutinated such blood was typed Rh positive.


The blood of the remaining 15% lacked the factor and was typed Rh negative.

Agglutination 85%


Rh positive RBCs

No agglutination 15%

Rh negative RBCs



LANDSTEINER



WIENER




LEVINE




Rh Blood Group System

Antigens of Rh systems: D d E e C c
more than 50 antigens



D antigen is on the RBC: **RhD positive**




D antigen lack on the RBC: **RhD negative**

d is only a symbol of lack the D antigen



RhD negative means there is no D antigen on the cell membrane



Basic of RhD grouping



Characteristics of IgM and IgG antibodies


IgM	IgG
	 HDN
e.g. ABO	e.g. Rh
Pentamer (big)	Monomer (small)
10 binding sites	2 binding sites
complete	incomplete
naturally occurring	Immuno antibody
NOT accross the placenta	ACCROSS the placenta
Reaction temperature +4 °C to room temperature	optimal reaction temperature 37 °C
Binds Complement	NOT binds Complement

Monoclonal IgM anti-D is used to in vitro tests

Reaction scheme of agglutination

- | | |
|------------------------------------|-----------------|
| 1. +++++ strong | positive |
| high, rough cloddy, hard frangible | |
| 2. +++ | variant |
| large, separated | |
| 3. ++ | |
| minor agglutination | negative |
| 4. + | |
| fine-grained | |
| 5. +/- | negative |
| very small | |
| 6. No agglutination, homogene | negative |



Anti-D + patient's RBC





INTERPRETATION OF RESULTS

anti-D


control




NO agglutination
Rh D negative



STRONG (++++) agglutination
Rh D positive



Weak agglutination (+ - +++)
Rh D varies




Negative control : anti-D control + patient's RBC
The control must be negative.
If any agglutination occurs the test is not valid.



What to do if the reaction is not ++++ positive?

1. Repeat the test
2. Send the sample to the Blood Bank
3. If the patient have to take transfusion he is RhD negative



Limitations of the Procedure:

Sources of errors

Antigen – antibody tests met a lot of requirements. Factors affecting the antigen-antibody reactions should be considered to establish the suitable reaction.

If the reaction conditions are not followed, false negative or false positive results can occur, which can lead to incorrect blood group determination.

- **Medium of reaction (ionic strength) – Serum**
- **Antigen - antibody ratio (50% suspension)**
- **Reaction temperature (+20 °C – room temperature)**
- **Reaction time 5 minutes**



Sympexis = rouleaux formation of RBCs

Physicochemical changes not real agglutination

The RBC's here have stacked together in long chains. This is known as "rouleaux formation" and it happens with increased serum proteins, particularly fibrinogen and globulins. Such long chains of RBC's sediment more readily. This is the mechanism for the sedimentation rate, which increases non-specifically with inflammation and increased "acute phase" serum proteins.



- Causes:**
- infections
 - multiple myeloma, cirrhosis
(an increase in the ratio of immunoglobulins to albumin)
 - inflammatory and connective tissue disorders
 - cancer
 - diabetes mellitus
 - an increase in the ratio of RBCs to plasma volume
(anemia, hypovolemia)
 - macromolecules, contrast medium

Acute phase proteins, particularly fibrinogen, interact with sialic acid on the surface of RBCs to facilitate the formation of rouleaux. Rouleaux formation is retarded by albumin proteins, in vitro by physiological saline.



Main causes of false positive reactions

- Rouleaux formation
marginal drying
- Little drops
- Late evaluation
- Contamination

Main causes of false negative reactions

- Early evaluation
- Inadequate
antigen-antibody ratio
- Expired reagents and test cells

The sympexis may be differentiate from real agglutination with dropping of phys. saline.

The sympexis dissolved but no agglutination.




TASKS

- RhD grouping with slide method – 2 tests


Two blood samples

1 empty tube respectively



- Bed side blood grouping

1 bed side card, one test




BLOOD SAMPLE LABELING

Details must include:


Patient's registered family name and given name (unabbreviated)

Date of birth or social security number (TAJ) (both preferred)




Name and code of institute requiring

(Signature or initials of the collector)
Date and time of collection



RhD typing – slide method

Labeling of slide



Data:

•Patient's name and code

•Reagent's name
Anti-D


•Negative control

Anti-D control
reagents + patient's
RBC



PREPARATION OF 50% RBC SUSPENSION


NATIVE BLOOD SAMPLE



RhD typing

Dropping of reagents

1 drop of



RhD typing

Dropping of RBC suspension

Kovács János 42.11.27. Date: 20013.11....

Anti-D control

Nagy Gizella 91.03.18.

Anti-D control

1 drpo of



- Patient's 50% RBC suspension to each reaction areas
- Mixing (wiping)
- Reaction time: 5 minutes
- Tilting
- Interpretation





INTERPRETATION OF RESULTS

anti-D

control




NO agglutination
Rh D negative



STRONG (++++) agglutination
Rh D positive



Weak agglutination (+ - +++)
Rh D varians



Negative control : anti-D control + patient's RBC


The control must be negative.

If any agglutination occurs the test is not valid.



Blood grouping on bedside card


ABO (forward) and RhD typing



Adhesive tapes



Stirring rods



Storage temperature: 2 – 25 °C

Lot number

Expiry date – important!



Expiry date

Lot number

Dried monoclonal specific antisera

Anti-A

Anti-B

Anti-D

Labelling of card


Name and codes

Date

Signature

Blood group result

Unit code
(RBC product grouping)




Autocontrol: physiological saline + patient's RBC



Test procedure

1 drop



Physiological saline
to the 4 reaction
fields respectively

Whole **blood** from
patient or blood unit

The blood may be native,
anticoagulated (**except
heparin**) venous or
capillary blood

Mixing until the reagent is completely dissolved

Spread material to be tested over the entire reaction field.

Repeat this procedure with a new or well cleaned mixing stick on the next reaction fields.



Blood being tested

Serum

	Anti-A	Anti-B	
Type AB (contains agglutinogens A and B)			AB
Type B (contains agglutinin B)			B
Type A (contains agglutinin A)			A
Type O (contains no agglutinogens)			O

AB

B

A

O



INTERPRETATION OF RESULTS

Anti-A	Anti-B	Anti-D	RESULT
positive	negative	negative	A RhD negative
positive	negative	positive	A RhD positive
negative	positive	positive	B RhD positive
negative	positive	negative	B RhD negative
positive	positive	positive	AB RhD positive
positive	positive	negative	AB RhD negative
negative	negative	negative	O RhD negative
negative	negative	positive	O RhD positive

Sources of errors:

Fals positive: rouleaux formation, drying out, cold reactive antibody

Fals negative: Blood drop too small or too large, expiry date has passed or card was stored improperly, red cell suspension (10%), hematocrit under 15 %




Bed side test

Interpretation of results

Anti-A	Anti-B	Anti-D	blood group
			A Rh pos
			A Rh neg
			B Rh pos
			B Rh neg
			AB Rh pos
			AB Rh neg
			O Rh pos
			O Rh neg



INTERPRETATION OF RESULTS




A RhD positive

Agglutination with anti-A and anti-D

Autocontrol should be negative



INTERPRETATION OF RESULTS




A RhD negative

Agglutination with anti-A only

Autocontrol should be negative



INTERPRETATION OF RESULTS




B RhD positive

Agglutination with anti-B and anti-D

Autocontrol should be negative



INTERPRETATION OF RESULTS



AB RhD negative

Agglutination with anti-A and anti-B

Autocontrol should be negative



INTERPRETATION OF RESULTS

Serafol® ABO+D

Ch.-B.Lot No. 080601 Expiry date: 31.12.01

Anti-A Anti-B Anti-D

Name (Name/ Nom) _____ ID _____
 Geb. Dat. (Date of Birth/ Date de Naissance) _____ Koms.- Nr. _____
 Datum (Date) _____ Blutgruppe (Blood Group/ Groupe Sanguin) _____
 Unterschrift (Signature) _____

Blut (Blood/Sang) **O RhD positive**

Produced by SIFIN GmbH D-13088 Berlin
 Distributed by Biotest AG D-63303 Dreieich

O RhD positive

Agglutination with anti-D only

Autocontrol should be negative

RECORD YOUR RESULT ON THE CARD



ORSZÁGOS VÉRELLÁTÓ SZOLGÁLAT

Serafol® ABO+D

Ch.-B.Lot No. 131107 Expiry date: 2008-11

Anti-A Anti-B Anti-D

Name (Name/ Nom) _____ ID _____
 Geb. Dat. (Date of Birth/ Date de Naissance) _____ Koms.- Nr. (Unit No./ No. Poche) _____
 Datum (Date) _____ Blutgruppe (Blood Group/ Groupe Sanguin) _____
 Unterschrift (Signature) _____

Blut (Blood/Sang) _____

Manufacturer: SIFIN GmbH D-13088 Berlin
 Distributor: Biotest AG D-63303 Dreieich

Serafol® ABO+D

Ch.-B.Lot No. 131107 Expiry date: 2008-11

Anti-A Anti-B Anti-D

Name (Name/ Nom) [redacted] ID 7890021
 Geb. Dat. (Date of Birth/ Date de Naissance) 20.02.69 Koms.- Nr. (Unit No./ No. Poche) _____
 Datum (Date) 01.08.2008 Blutgruppe (Blood Group/ Groupe Sanguin) AB Rh(+)
 Unterschrift (Signature) _____

Blut (Blood/Sang) _____

Manufacturer: SIFIN GmbH D-13088 Berlin
 Distributor: Biotest AG D-63303 Dreieich

Serafol® ABO+D

Ch.-B.Lot No. 131107 Expiry date: 2008-11

Anti-A Anti-B Anti-D

Name (Name/ Nom) _____ ID _____
 Geb. Dat. (Date of Birth/ Date de Naissance) _____ Koms.- Nr. (Unit No./ No. Poche) _____
 Datum (Date) _____ Blutgruppe (Blood Group/ Groupe Sanguin) _____
 Unterschrift (Signature) _____

Blut (Blood/Sang) _____

Manufacturer: SIFIN GmbH D-13088 Berlin
 Distributor: Biotest AG D-63303 Dreieich

Serafol® ABO+D

Ch.-B.Lot No. 131107 Expiry date: 2008-11

Anti-A Anti-B Anti-D

Name (Name/ Nom) [redacted] ID 7890053
 Geb. Dat. (Date of Birth/ Date de Naissance) 10.02.70 Koms.- Nr. (Unit No./ No. Poche) M101151321007 21 X
 Datum (Date) 01.08.2008 Blutgruppe (Blood Group/ Groupe Sanguin) AB Rh(+)
 Unterschrift (Signature) _____

Blut (Blood/Sang) _____


Manufacturer: SIFIN GmbH D-13088 Berlin
 Distributor: Biotest AG D-63303 Dreieich



Antibody detection

Immune antibodies are IgG origin

Transfusions, pregnancy, transplantation




Clinical importance of immune antibodies:

May cause agglutination or hemolysis in vivo – maybe fatal

Hemolytic disease of the newborn

Type of antibodies




ORSZÁGOS VÉRELLÁTÓ SZOLGÁLAT



Coombs reaction

Anti-human globulin (Coombs) antibody are produced by immunizing non-human species with human serum.


Direct Coombs test / Direct antiglobulin test



Animal anti-human antibodies will also bind to human antibodies, commonly IgG or IgM that may be fixed onto antigens on the surface of red blood cells and in the appropriate test conditions this can lead to agglutination of RBCs.

The phenomenon of agglutination of RBCs is important here, because the resulting clumping of RBCs can be visualised; **when clumping is seen the test is positive and when clumping is not seen the test is negative.**

Indirect Coombs test / Indirect antiglobulin test




The DAT is used to detect IgG or C3 bound to the surface of the red cell.

The IAT is used to detect free red cell antibodies in patient serum.

Antibody tests

Test methods for antibody detection

Tube test, microplate test, **column agglutination**



Gel separates agglutinates based on size and by binding to IgG-coated RBCs
 a) Strong positive gives agglutinates at the TOP of the gel (left side of above)
 b) Complete negative gives RBCs at the BOTTOM of the gel (right side)

Test procedure


Patient's serum (**unknown antibody**) + test cells (**known antigens**)

Reaction temperature + **37 oC**



Importance of antibody identification

- HDN
- Transfusion reaction
- Compatibility



Types of the test

Antibody screening




ID-DiaCell I-II-III/Antigen Table (example)

Rh-hr Donor	Rh						Kell						Duffy		Kidd		Lewis		P	MNS				Luth.		Xg		
	D	C	E	c	e	C ^w	K	k	Kp ^a	Kp ^b	Js ^a	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a	F/M	
I C ^w CD.ee R ₁ ^w R ₁ 754456	+	+	0	0	+	+	0	+	0	+	0	+	+	+	0	+	+	0	+	+	0	+	0	+	0	+	0	M
II ccD.EE R ₂ R ₂ 047144	+	0	+	+	0	0	0	+	+	+	0	+	+	0	+	+	0	+	0	+	+	+	+	0	+	+	+	M
III ccddee rr 101469	0	0	0	+	+	0	+	+	0	+	0	+	0	+	+	0	0	+	+	0	+	0	+	+	+	+	+	F



Antibody identification



Rh-hr	Donor	Rh						Kell						Duffy		Kidd		Lewis		P	MNS				Luth.		Xg		
		D	C	E	c	e	C ^w	K	k	Kp ^a	Kp ^b	Js ^a	Js ^b	Fy ^a	Fy ^b	Jk ^a	Jk ^b	Le ^a	Le ^b	P ₁	M	N	S	s	Lu ^a	Lu ^b	Xg ^a	F/M	
1	C ^w CD.ee R ₁ W ₁ 677783	+	+	0	0	+	+	0	+	0	+	+	0	+	0	0	+	+	+	+	0	+	0	0	+	+			
2	CCD.ee R ₁ R ₁ 113683	+	+	0	0	+	0	+	+	0	+	0	+	0	+	+	0	0	+	0	+	0	+	0	0	+	0		
3	ccD.EE R ₂ R ₂ 422278	+	0	+	+	0	0	0	+	0	+	0	+	+	+	0	+	0	+	+	+	+	+	0	0	+	+		
4	Ccddee r'r 293832	0	+	0	+	+	0	0	+	0	+	0	+	+	0	0	+	+	+	+	+	+	+	0	0	+	+		
5	ccddEe r''r 307849	0	0	+	+	+	0	0	+	0	+	0	+	+	+	+	0	+	+	+	+	+	0	+	+	+	0		
6	ccddeEe rr 308478	0	0	0	+	+	0	+	+	0	+	0	+	0	+	0	0	+	0	+	0	+	0	+	0	+	nt		
7	ccddeEe rr 439656	0	0	0	+	+	0	0	+	0	+	0	+	0	+	0	+	+	0	+	+	0	+	+	0	+	+		
8	ccD.ee Ror 032656	+	0	0	+	+	0	0	+	0	+	0	+	0	0	+	0	0	+	0	+	0	+	0	+	0	+		
9	ccddeEe rr 341656	0	0	0	+	+	0	0	+	0	+	0	+	0	0	+	0	+	+	0	0	+	0	+	0	+	0		
10	ccddeEe rr 454253	0	0	0	+	+	0	0	+	0	+	0	+	0	0	+	0	+	+	0	0	+	0	+	0	+	+		
11	ccddeEe rr 169276	0	0	0	+	+	0	0	+	+	0	+	+	0	+	0	+	+	+	0	0	+	0	+	0	+	+		



Examples: Anti-D



N°	GENOTYPE	Rhesus						Kell				Duffy		Kidd		Lewis		MNS				P	Lutheran		Colton *			
		D RH1	C RH2	E RH3	c RH4	e RH5	Cw RH8	K KEL1	k KEL2	Kp ^a KEL3	Kp ^b KEL4	Fy ^a FY1	Fy ^b FY2	Jk ^a JK1	Jk ^b JK2	Le ^a LE1	Le ^b LE2	M MNS1	N MNS2	S MNS3	s MNS4	P1 P1	Lu ^a LU1	Lu ^b LU2	Coa	Cob		
IIP	R1wR1 : RH(1,2,-3,-4,5,8)	+	+	0	0	+	+	0	+	0	+	0	+	0	+	0	+	+	0	0	+	0	0	+	0	+	nt	0
IIIP	R2R2 : RH(1,-2,3,4,-5,-8)	+	0	+	+	0	0	0	+	0	+	+	+	+	+	0	+	+	0	+	+	+	+	+	0	nt	0	
IIIP	rr : RH(-1,-2,-3,4,5,-8)	0	0	0	+	+	0	+	+	0	+	0	+	+	0	0	+	0	+	0	+	0	+	0	+	nt	0	


N°	GENOTYPE	Rhesus						Kell				Duffy		Kidd		Lewis		MNS				P	Lutheran						
		D RH1	C RH2	E RH3	c RH4	e RH5	Cw RH8	K KEL1	k KEL2	Kp ^a KEL3	Kp ^b KEL4	Fy ^a FY1	Fy ^b FY2	Jk ^a JK1	Jk ^b JK2	Le ^a LE1	Le ^b LE2	M MNS1	N MNS2	S MNS3	s MNS4	P1 P1	Lu ^a LU1	Lu ^b LU2					
1/1P	R1wR1 : RH(1,2,-3,-4,5,8)	+	+	0	0	+	+	0	+	0	+	+	+	+	0	0	+	+	0	+	0	+	0	+	0	+			
2/2P	R1wR1 : RH(1,2,-3,-4,5,8)	+	+	0	0	+	+	0	+	0	+	+	0	0	+	0	+	+	+	0	+	+	0	+	+	0	+		
3/3P	R1R2 : RH(1,2,3,-4,5,-8)	+	+	+	0	+	0	0	+	0	+	+	+	+	0	0	+	+	+	0	+	+	0	+	+	0	+		
4/4P	R2R2 : RH(1,-2,3,4,-5,-8)	+	0	+	+	0	0	0	+	0	+	0	+	0	+	+	0	+	+	0	+	+	0	+	+	0	+		
5/5P	R0r : RH(1,-2,-3,4,5,-8)	+	0	0	+	+	0	+	+	0	+	+	0	+	+	0	+	+	0	+	+	0	+	+	0	+			
6/6P	r'r : RH(-1,2,-3,4,5,-8)	0	+	0	+	+	0	0	+	0	+	+	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+		
7/7P	r''r : RH(-1,-2,3,4,5,-8)	0	0	+	+	+	0	0	+	0	+	0	+	+	0	0	+	+	+	+	+	+	+	+	+	+			
8/8P	rr : RH(-1,-2,-3,4,5,-8)	0	0	0	+	+	0	+	+	0	+	+	+	0	+	0	+	0	+	0	+	0	+	0	+	0	+		
9/9P	rr : RH(-1,-2,-3,4,5,-8)	0	0	0	+	+	0	0	+	+	0	+	+	0	+	0	+	+	+	+	+	+	+	+	0	0	+		
10/10P	rr : RH(-1,-2,-3,4,5,-8)	0	0	0	+	+	0	0	+	0	+	0	+	0	+	0	+	+	+	+	+	+	+	0	0	+			



Examples:

Anti-D or
Anti-Jkb or ?

Anti-M




N°	Antigen table GENOTYPE	Rhesus						Kell				Duffy		Kidd		Lewis		MNS				P	Lutheran		Colton *		
		D RH1	C RH2	E RH3	c RH4	e RH5	Cw RH8	K KEL1	k KEL2	Kp ^a KEL3	Kp ^b KEL4	Fy ^a FY1	Fy ^b FY2	Jk ^a JK1	Jk ^b JK2	Le ^a LE1	Le ^b LE2	M MNS1	N MNS2	S MNS3	s MNS4	P1 P1	Lu ^a LU1	Lu ^b LU2	Coa Coa	Cob Cob	
IP	R1wR1 : RH(1,2,-3,-4,5,8)	+	+	0	0	+	+	0	+	0	+	0	+	0	+	+	0	0	+	0	0	+	0	0	+	nt	0
IIP	R2R2 : RH(1,-2,3,4,-5,-8)	+	0	+	+	0	0	0	+	0	+	+	+	+	+	0	+	+	0	+	+	+	+	0	nt	0	
IIIP	rr : RH(-1,-2,-3,4,5,-8)	0	0	0	+	+	0	+	+	0	+	0	+	0	0	+	0	+	0	+	+	+	0	+	nt	0	
PATIENT																											

ORSZAGOS VERELLATÓ SZOLGALAT



Examples:

Anti-M



N°	Panel table GENOTYPE	Rhesus						Kell				Duffy		Kidd		Lewis		MNS				P	Lutheran		
		D RH1	C RH2	E RH3	c RH4	e RH5	Cw RH8	K KEL1	k KEL2	Kp ^a KEL3	Kp ^b KEL4	Fy ^a FY1	Fy ^b FY2	Jk ^a JK1	Jk ^b JK2	Le ^a LE1	Le ^b LE2	M MNS1	N MNS2	S MNS3	s MNS4	P1 P1	Lu ^a LU1	Lu ^b LU2	
1/1P	R1wR1 : RH(1,2,-3,-4,5,8)	+	+	0	0	+	+	0	+	0	+	+	+	+	0	0	+	+	0	+	0	+	+	0	+
2/2P	R1wR1 : RH(1,2,-3,-4,5,8)	+	+	0	0	+	+	+	+	0	+	+	0	0	+	0	+	+	+	0	+	+	0	+	
3/3P	R1R2 : RH(1,2,3,-4,5,-8)	+	+	+	0	+	0	0	+	0	+	+	+	+	0	0	+	+	+	0	+	+	0	+	
4/4P	R2R2 : RH(1,-2,3,4,-5,-8)	+	0	+	+	0	0	0	+	0	+	0	+	+	0	+	+	+	0	+	+	+	0	+	
5/5P	R0r : RH(1,-2,-3,4,5,-8)	+	0	0	+	+	0	+	+	0	+	+	0	+	+	0	+	+	+	0	+	+	0	+	
6/6P	r'r : RH(-1,-2,-3,4,5,-8)	0	+	0	+	+	0	0	+	0	+	+	+	0	0	+	0	+	0	+	0	+	0	+	
7/7P	r'r : RH(-1,-2,3,4,5,-8)	0	0	+	+	+	0	0	+	0	+	+	+	0	0	+	+	+	+	0	+	+	0	+	
8/8P	rr : RH(-1,-2,-3,4,5,-8)	0	0	0	+	+	0	+	+	0	+	+	+	+	0	0	+	0	+	0	+	+	0	+	
9/9P	rr : RH(-1,-2,-3,4,5,-8)	0	0	0	+	+	0	0	+	+	+	0	+	+	0	0	+	+	+	+	+	0	0	+	
10/10P	rr : RH(-1,-2,-3,4,5,-8)	0	0	0	+	+	0	0	+	0	+	+	0	+	0	0	+	+	+	+	0	+	0	+	

ORSZAGOS VERELLATÓ SZOLGALAT

